

SAMPLING KITS



WATER SAMPLING KIT

The Water Sampling Kit is designed for the aseptic collection, transport and analysis of chlorinated and non-chlorinated water samples. The kit includes a non-toxic neutralizing tablet of sodium thiosulfate to enable detection of potential pathogens. Approved by the UPA for drinking water analysis, it is an economical and efficient replacement for rigid sampling containers.

- Analysis of water from different sources
- Sodium thiosulfate tablet
- Volume markings

Indicated volume = Functional volume - Corresponds to the fill line.



ECOLO

Standard closure PRINTING - WRITING AREA WITH GRADUATION

SKU	Solution	Dimension		Volume		Thickness		Packaging
		in x in	mm x mm	ml	oz	mils	µm	
KWS-21100	10 mg	3.5 x 7	89 x 178	100	4	3	76	100 (1 x 100)
KWS-22200	30 mg	5 x 9	127 x 229	300	12	3	76	100 (1 x 100)

Heavy sample closure PRINTING - WRITING AREA WITH GRADUATION

SKU	Solution	Dimension		Volume		Thickness		Packaging
		in x in	mm x mm	ml	oz	mils	µm	
KWS-33100	10 mg	5 x 12	127 x 305	500	20	4	101	100 (1 x 100)
KWS-34200	30 mg	7 x 12	178 x 305	1000	35	4	101	100 (1 x 100)

Ecolo - Standard closure PRINTING - WRITING AREA WITH GRADUATION

SKU	Solution	Dimension		Volume		Thickness		Packaging
		in x in	mm x mm	ml	oz	mils	µm	
KWS-21100E	10 mg	3.5 x 7	89 x 178	100	4	3	76	100 (1 x 100)
KWS-22200E	30 mg	5 x 9	127 x 229	300	12	3	76	100 (1 x 100)

Ecolo - Heavy sample closure PRINTING - WRITING AREA WITH GRADUATION

SKU	Solution	Dimension		Volume		Thickness		Packaging
		in x in	mm x mm	ml	oz	mils	µm	
KWS-33100E	10 mg	5 x 12	127 x 305	500	20	4	101	100 (1 x 100)
KWS-34200E	30 mg	7 x 12	178 x 305	1000	35	4	101	100 (1 x 100)



Inspired by your challenges